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
1. Adult female Swiss mice were sacrificed by cervical dislocation. Lumbar dorsal root ganglia were collected in sterile conditions.
2. Ganglia were successively treated by two incubations with collagenase A (1 mg/ml, Roche Diagnostic) for 45 minutes each (37 °C) and then with trypsin-EDTA (0.25 %, Sigma) for 30 minutes.
3. They were mechanically dissociated by passing 8 to 10 times through the tip of a fire-polished Pasteur pipette in neurobasal (Life Technologies) culture medium supplemented with 10% foetal bovine serum (FBS).
4. Isolated cells were collected by centrifugation and suspended in neurobasal culture medium not supplemented with FBS.
5. Dissociated neurons were counted with a Malassez cell and plated on D,L-polyornithine (0.5 mg/ml)-laminin (5 µg/ml)-coated glass coverslips at a density of 2500 neurons per coverslip, in 4-wells culture-treated polystyrene plates with neurobasal culture medium not supplemented with FBS, 2 mM glutamine, penicillin/streptomycin (20 U/ml, 0.2 mg/ml).
6. Plates were incubated in an incubator with a humidified 95 % air/ 5 % CO<sub>2</sub> atmosphere.

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### 5a : Protocoles

Mouse peritoneal macrophages in RPMI 1640 medium supplemented with 10 % fetal bovine serum (FBS) were seeded in one dish with PNIPAM Surface and two traditional cultureware dishes at  $2.4 \times 10^5$  cells/cm<sup>2</sup>.

The cells were incubated at 37 °C in a humidified atmosphere of 5 % CO<sub>2</sub> in air.


Cells were then cultured for 2 days in RPMI 1640 medium supplemented with 10 % FBS and harvested using one of the following procedures:

| Protocol 1 : temperature reduction on PNIPAM surface                                                                                                                                                                               | Protocol 2 : trypsinization                                                                                                                                                                                                                                                     | Protocol 3 : scraping*                                                                                                                                                                                                     |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1. Cells were washed with Ca <sup>2+</sup> - and Mg <sup>2+</sup> -free PBS<br>2. 4.0 mL RPMI 1640 medium supplemented with 10 % FBS was added and the dish was incubated at 20 °C for 30 min.<br>3. Detached cells were harvested | 1. Cells were washed with Ca <sup>2+</sup> - and Mg <sup>2+</sup> -free PBS<br>2. 1.0 ml of 0.25% trypsin/EDTA was added, and the dish was incubated at 37 °C for 5 min.<br>3. 3.0 ml RPMI 1640 medium supplemented with 10 % FBS was added<br>4. Detached cells were harvested | 1. Cells were washed with Ca <sup>2+</sup> - and Mg <sup>2+</sup> -free PBS<br>2. 4.0 ml of 2.5 mM EDTA/PBS was added and the dish was incubated on ice for 20 min<br>3. The cells were detached by scraping and harvested |

The harvested cells were counted and the recovery ratio ((detached cells / total cells) x 100) was calculated.

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
### 6b : Résultats

**Grey** : anti-CD140a

**White** : isotype control for non-specific binding (background)

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Copyright © 2026 FormaV. Tous droits réservés.